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(IE, MT, UK) Safety Instructions

Please take a little time and read through these safety instruction carefully. We would like to provide you with important information on the safe use of PILZ® laboratory heating mantles. This applies above all to your own protection in the event of an unexpected reaction in the course of heating liquids, overflowing or a malfunction. The type plate informs you about the technical data and connected loads of the device.

Store these operating instructions within easy reach near the place of use. Please also observe the guidelines for ensuring occupational safety as well as those issued by the employers' liability insurance associations (i.e. TRGS526 - Technical Regulations for Hazardous Substances Laboratories). Include these notes as part of work instructions in your quality management manual.

Also treat them as operating instructions according to TRGS526 section 4.2.

Warning labels on the product must not be removed! The PILZ® laboratory heating mantle is electrical equip-

ment.

Maintenance is the responsibility of a qualified electrician, or a person trained in electrical engineering.

In addition, PILZ® laboratory heating mantles are equipped with a residual current circuit breaker in the connection plug. If an electrician removes the residual current circuit breaker, the electrician must ensure that the mains supply line is protected by a residual current circuit breaker.

If an electrical fault occurs, it safely disconnects the laboratory heating hood from the operating voltage and prevents personal injury.

Before commissioning and at regular intervals, the residual current circuit breaker must be checked for proper functioning.

Moisture protection of the laboratory heating mantle

When using the PILZ® heating jacket, the use of a residual current circuit breaker in accordance with the standards and regulations is mandatory. If electrically conductive liquid penetrates the heating element and it becomes soaked, the residual current circuit breaker integrated in the power supply line or connection plug switches reliably and permanently voltage-free (switching threshold: residual current >30 mA [DIN EN 61008-1]). If this occurs, the laboratory heating mantle must be disconnected from the mains immediately by pulling the connection plug out of the protective contact socket.

**The PILZ® Protect Line laboratory heating mantle is additionally equipped with a heating element with metallic sheathing as a protective conductor. The design guarantees the operating mode of protection class I for the PILZ® laboratory heating mantle and triggers the residual current circuit breaker more quickly.

The PILZ® laboratory heating mantle must not be operated after the protective device has been triggered!

Carry out a function and safety check before operation. The soaked/electrically faulty PILZ® laboratory heating mantle must be marked as „blocked“ and sent to the blocked storage. The PILZ® laboratory heating mantle must be stored in a drying cabinet, i.e. for drying. The laboratory heating mantle must be checked for dryness before being put back into operation. It must not be put into operation until it has been approved by the responsible electrician or electrically instructed person.

The PILZ® laboratory heating mantle must also be taken out of operation if,

- Protect Line Series; the protective conductor (braiding) of the heating element is damaged.
- Wire connection lines show visible damage (i.e. cuts, abrasion, brittle and cracked sheathing).
- Loose plugs and cable connections can be detected.
- Visible heating cable wire in the interior (mechanical damage).
- brittle or broken glass fabric (overheating of the glass fabric).

A qualified electrician or person trained in electrical engineering decides on possible repair or disposal.

Prevent overheating!

To always ensure safe operation of the PILZ® device heating mantle, overheating must be avoided. This can be achieved, for example, by using a temperature controller. Also, in the case of devices with heating stage switches, only those stages should operate in heating mode for which a liquid level can be detected in the glass bulb. In the case of heating mantles with power switching, the heating due to the reaction of the liquid and the filling level must be observed to avoid overheating!

For products with multiple heating zones, as the liquid

level drops, the heating circuits that are not necessary must be switched off.

PILZ® laboratory heating mantles must not be operated with empty, encrusted or completely without glass bulbs.

Do not touch the PILZ® heating mantle when it is switched on or ready for operation.

Do not touch the heating surface with metal or electrically conductive objects.

Do not use metal vessels, vessels with metal parts or electrically conductive vessels in the heating surface area. Do not insert needles or similar into the PILZ® heating mantle, as this may damage heating coils or other electrical conductors and/or their protective devices. Electrical shocks with serious injuries may result.

The PILZ® heating mantle contains a heating element in the interior. This heating element must not be touched when the PILZ® heated mantle is switched on or still hot. The upper edge of the PILZ® heated mantle and the fastening hooks can reach high temperatures. There is a risk of injury from burns.

PILZ® laboratory heating mantles are not approved for use in potentially explosive atmospheres!

It must not be used for heating processes where danger may arise from explosive media or explosive gas-air mixtures.

Note that fire can also occur with flammable liquids in the event of flask breakage, overfilling and spillage.

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Winkler AG